



Performance
Made Compact

The next generation of Dynamic Support Technology™

The F&P Evora™ Full mask is uniquely designed to optimise performance throughout the night.

Dynamic Support Technology combines Stability Wings with the Floating Seal to allow for freedom of movement while keeping the mask in place.

Stability



rated the Evora Full mask as stable or very stable¹

Performance



rated the Evora Full seal performance as very good or good²



Compact Full Face Mask



Headgear with VentiCool™ Technology

Breathable headgear designed to allow for heat and moisture to escape

Removable Frame Clip

Enables the headgear settings to stay in place during fitting and cleaning

Exhaust Holes

Designed for minimal noise

Considering patient preference

Up to 60% of PAP therapy users report issues such as poor mask fit, air leak, redness or skin irritation, and claustrophobia.³⁻⁹

When choosing a new mask, consider the mask design and patient preference.



never or rarely experienced unintentional air leak into their eyes while using the Evora Full mask¹⁰



never or rarely felt claustrophobic while using the Evora Full mask¹¹

F&P Evora Full

Sizing options



XS


Extra Small

S-M

Small to Medium

L

Large



Std.

Standard

XL

Extra Large

Fitting summary



1. Adjust the top strap



2. Adjust the bottom strap



3. Check the space

For further instruction, please refer to the Fitting Your Mask section of the Evora Full Use and Care Guide.

Product Codes

Code	Description
EVFIXU	Evora Full Mask - Extra Small U Model
EVFIMU	Evora Full Mask - Small to Medium U Model
EVFILU	Evora Full Mask - Large U Model
EVFIXXU	Evora Full Mask - Double Seal Extra Small U Model
EVFIMMU	Evora Full Mask - Double Seal Small to Medium U Model
EVFILLU	Evora Full Mask - Double Seal Large U Model
EVFIXMLU	Evora Full Mask Fitpack - Extra Small, Small to Medium, Large U Model

Spare Parts

Code	Description
400EVF111	Evora Full Seal - Extra Small Spare
400EVF112	Evora Full Seal - Small to Medium Spare
400EVF113	Evora Full Seal - Large Spare
400EVF121	Evora Full Headgear Spare - Standard
400EVF122	Evora Full Headgear Spare - Extra Large
400EVF151	Evora Full Tube Spare
400EVF161	Evora Full Frame Spare
400FPH341	F&P Swivel White

F&P Evora Full is available on the F&P myMask™ App

Instructs patients on mask use through step-by-step fitting, fine-tuning and cleaning videos. Find out more. Scan the QR code, or download directly.



1. 40 out of 44 participants rated the Evora Full stability as either very stable or stable. Internal validation study conducted with 44 participants in New Zealand 2020. 2. 39 out of 43 participants rated the overall seal performance of the Evora Full mask as very good or good. Internal validation study conducted with 44 participants in New Zealand 2020. 3. Ewles et al 2018. doi:10.1183/13993003.congress-2018.pa2372. 4. Patout et al 2018. DOI: 10.1183/13993003.congress-2018.PA2378. 5. Andreu et al 2012. doi:10.1183/09031936.00013311. 6. Brown et al 2013. doi:10.1016/j.jstrokecerebrovasdis.2011.06.010. 7. Nosedo et al 2009. doi:10.1016/j.sleep.2008.08.010. 8. Masa et al 2004. doi:10.1164/rccm.200312-1787OC. 9. Nussbaumer et al 2006. doi:10.1378/chest.129.3.638. 10. 41 out of 44 participants never or rarely experienced unintentional mask leak into the eyes when sleeping with the Evora Full mask. Internal validation study conducted with 44 participants in New Zealand 2020. 11. 43 out of 44 participants never or rarely felt claustrophobic while wearing the Evora Full mask. Internal validation study conducted with 44 participants in New Zealand 2020.

iPhone is a trademark of Apple Inc., registered in the U.S. and other countries.
F&P, Evora, myMask, Dynamic Support Technology and VentiCool are trademarks of Fisher & Paykel Healthcare Limited.
For patent information, please see www.fphcare.com/ip. REF 625323 REV A © 2021-11 Fisher & Paykel Healthcare Limited.