

F&P Solo mask is intended to be used for single patient reuse in the home and multi-patient re-use in a hospital or other clinical setting.



Characteristics of the masks	F&P Solo Nasal	F&P Solo Pillows
Type of mask	Under-the-nose nasal	Under-the-nose with prongs
Contraindications	None	None
Storage temperature	-20 °C to 50 °C (-4 °F to 122 °F)	-20 °C to 50 °C (-4 °F to 122 °F)
Operating temperature	5 to 40 °C (41 to 104 °F).	5 to 40 °C (41 to 104 °F).
Operating pressure	4 to 20 cmH ₂ O	4 to 20 cmH ₂ O
Flow	15 to 49 L/min	15 to 49 L/min
Sound power level	The sound power level of the mask is 31.5 dBA, with uncertainty ± 2.5 dBA.	The sound power level of the mask is 31.5 dBA, with uncertainty ± 2.5 dBA.
Sound pressure level	The sound pressure level of the mask is 23.6 dBA, with uncertainty ± 2.5 dBA	The sound pressure level of the mask is 23.6 dBA, with uncertainty ± 2.5 dBA
Interface connections	ISO 5356-1 Conical Connectors	ISO 5356-1 Conical Connectors
Resistance to flow through mask	Pressure drop through the nasal mask at 50L/min: 0.8 ± 0.5 cmH ₂ O. Pressure drop through the nasal mask at 100 L/min: 2.9 ± 1.0 cmH ₂ O.	Pressure drop through the pillows mask at 50 L/min: 1.6 ± 1.1 cmH ₂ O. Pressure drop through the pillows mask at 100 L/min: 6.0 ± 3.1 cmH ₂ O.
Quantities per box	10	10
Shelf life	2 Years	2 Years

Materials	Components	Materials
	Frame	Nylon
	Tube	Polycarbonate (PC)/Polyurethane
	Swivel	Nylon
	Seal	Polycarbonate (PC)/ Silicone
	Headgear	Poyether Block Amide (PEBA)/Nylon/ Silicone/Polycarbonate (PC)/Carbon Reinforced Nylon

Biocompatibility	F&P Solo Nasal	F&P Solo Pillows
	Meets standards: ISO 10993-1:2020	Meets standards: ISO 10993-1:2020
Packaging materials	Outer box: cardboard In the package: Sugar cane pulp tray	Outer box: cardboard In the package: Sugar cane pulp tray
Packaging dimensions (mm)	Length: 170, Height: 287, Width: 50	Length: 170, Height: 287, Width: 50
Disposal	The F&P Solo mask does not contain any hazardous materials and may be disposed of with general waste	The F&P Solo mask does not contain any hazardous materials and may be disposed of with general waste

Regulatory	F&P Solo Nasal	F&P Solo Pillows
Classification	USA - Class II; EU - Class IIa; Australia - Class IIa	USA - Class II; EU - Class IIa; Australia - Class IIa
Country of origin	New Zealand	New Zealand
GMDN code	57815	57815
CE marking	CE0123	CE0123
EU notified Body	TÜV 0123 - Product Service GMBH	TÜV 0123 - Product Service GMBH

Product Order Codes			
Code - Solo Nasal	Description	Code - Solo Pillow	Description
SLN1SU, SLN1SSU*	F&P Solo Nasal Mask Small U Model	SLP1SU, SLP1SSU*	F&P Solo Pillows Mask Small U Model
SLN1MU, SLN1MMU*	F&P Solo Nasal Mask Medium U Model	SLP1MU, SLP1MMU*	F&P Solo Pillows Mask Medium U Model
SLN1LU, SLN1LLU*	F&P Solo Nasal Mask Large U Model	SLP1LU, SLP1LLU*	F&P Solo Pillows Mask Large U Model
SLN1WU	F&P Solo Nasal Mask Wide U Model	SLP1SMLU	F&P Solo Pillows Mask FitPack U Model
SLN1SMLU	F&P Solo Nasal Mask FitPack U Model		

Spare Part Order Codes			
Code - Solo Nasal	Description	Code - Solo Pillow	Description
400SLN111	F&P Solo Nasal Cushion Small	400SLP111	F&P Solo Pillows Cushion Small
400SLN112	F&P Solo Nasal Cushion Medium	400SLP112	F&P Solo Pillows Cushion Medium
400SLN113	F&P Solo Nasal Cushion Large	400SLP113	F&P Solo Pillows Cushion Large
400SLN114	F&P Solo Nasal Cushion Wide		
400SL121	F&P Solo Headgear Spare		
400SL141	F&P Solo Tube, Frame & Swivel Spare		

GTIN		
Nasal Code	Each	Box
SLN1SU	09420012468222	09420012471352
SLN1MU	09420012468239	09420012471253
SLN1LU	09420012468246	09420012471215
SLN1SSU	09420012468277	09420012471345
SLN1MMU	09420012468284	09420012471239
SLN1LLU	09420012468291	09420012471192
SLN1WU	09420012468253	09420012471390
SLN1SMLU	09420012468260	09420012471314
400SLN111	09420012468550	09420012470775
400SLN112	09420012468567	09420012470782
400SLN113	09420012468574	09420012470799
400SLN114	09420012468581	09420012470805
400SL121	09420012468772	09420012470751
400SL141	09420012468789	09420012470768
Pillow Code	Each	Box
SLP1SU	09420012468307	09420012471581
SLP1SSU	09420012468345	09420012471574
SLP1MU	09420012468314	09420012471499
SLP1MMU	09420012468352	09420012471475
SLP1LU	09420012468321	09420012471444
SLP1LLU	09420012468369	09420012471420
SLP1SMLU	09420012468338	09420012471543
400SLP111	09420012468598	09420012471093
400SLP112	09420012468604	09420012471109
400SLP113	09420012468611	09420012471116
400SLP114	09420012468802	09420012471123